



2 Witney Way
Boldon Business Park
Tyne and Wear
NE35 9PE
England

Bipolar cables and Connectors

The following products are now available from Altomed Limited.



To place an order contact Customer Services on
Tel: (0) 191 519 0111 Fax: (0) 191 519 0283

Bipolar Machine End

British Standard (BS) #1
Fits old Downs, Siemens, Erbe



Half Round
Fits old Eschmann, G.U.



Two Plug
Twin 4mm pin, fits most modern machines



Codman Fitting
Butterfly plug with two pins



Bipolar Instrument End

British Standard (BS) #2
Accepts block end Eschman, Aesculap, etc.



2 Pin
Accepts twin round pins, Codman etc.



Two Flat Pin
Accepts block end Eschmann, Aesculap, etc.



Storz Fitting
Small silver plug with central pin



Altomed Standard Cables

A5080: 2 Plugs to 2 Pin
A5081: 2 Pin to 1/2 Round
A5082: BS#1 to 2 Pin
A5083: MIRA to 2 Pin
A5084: BS#1 to BS#2

Altomed Standard Cables

A5085: 1/2 Round to BS#2
A5086: 2 Plugs to BS#2
A5090: 2 Plugs to 2 Flat Pin
A5091: 1/2 Round to 2 Flat Pin
A5092: BS#1 to 2 Flat Pin

A standard cable is 3 meters (10') long. Any length or connection variation can be supplied. All cables are fully Autoclavable.