



2 Witney Way  
Boldon Business Park  
Tyne and Wear  
NE35 9PE  
England

# The Original Lester Jones and modified Implantable Tubes

The following products are now available from Altomed Limited.




To place an order contact Customer Services on  
Tel: (0) 191 519 0111 Fax: (0) 191 519 0283

Due to the multitude of different sizes and styles available, it is not practical to quote every single code number here therefore when ordering a tube, please quote the necessary code number by using the format outlined below. The Purchas Code is made by quoting first the Style (see below in red) then the Flange Width and the Length. If a Suture Hole is required add H to the end of the code, if it needs to be Angled add the letter A, and for a Frosted Finish add F. For example:

Examples of Products Required	Style	Flange Size	Tube Length	With Suture Hole	Angled	Frosted	Purchase Code
Lester Jones 3.5mm flange 9mm long with hole and angled	LJ	35	09	H	A		LJ3509HA
Gladstone Putterman 3.0mm flange 20mm long with hole and straight	GP	30	20	H			GP3020H
Naugle Tube 4.0mm flange 25mm long no hole and angled	NG	40	25		A		NG4025A
Callahan Cox Tube 3.0mm flange 28mm long no hole and straight	CC	30	28				CC3028
Lester Jones 3.5mm flange 15mm long, hole, straight and frosted finish	LJ	35	15	H		F	LJ3515HF

Please note: All Glass Implantable Tubes are supplied non-sterile and require sterilising in the autoclave before use.


**LJ** - Lester Jones Tube, Straight and Angled  
(2.2mm Outside Diameter Tube)



Jones Gold Dilator




A11908S – Polyethylene Tube  
Supplied Sterile and cut to length by surgeon



**GP** - Gladstone Putterman Tube



**BY** – Baylis 60° Curved Tube



**DC** – Double Collar Tube



**CC** – Callahan Cox Tube  
(3.0mm Outside Diameter Tube)



**NG** – Naugle Tube



## Related Devices

- A11937 – Lester Jones Gold Dilator No. I: 1.7mm and 1.9mm arms
- A11937A – Lester Jones Gold Dilator No. II: 2.2mm and 2.4mm arms
- A11937B – Lester Jones Gold Dilator No. III: 2.7mm and 2.9mm arms
- A11934 – Blue anodized aluminium identification scale for tubes

Ask for copies of the **Parkin Lester Jones Tube Introducer (Promo01)** and **Cox System (Promo03)** Promotional Literature to see the instruments designed to help with the implantation of the tubes.

Contact Altomed for pricing and availability or any further details.