





New CE marked Single Use Alger Brush Rust Ring Remover Burrs.

The Burrs are sold Sterile and ready for use



NEW →



Alger Brush Rust Ring Remover

During metal working procedures such as welding, sanding or grinding, it is possible that a spot of rust can be ejected from the machine into the patient's cornea. If this occurs it is essential that all the rust is removed as any remnants can cause irritation for days or even weeks. The Alger Brush has been specially designed to remove rust spots from a patient's cornea.

Special features of the Alger Brush:

- ✓ Motor will stop before injuring the cornea if too much pressure is applied
- ✓ The gentle slow rotation rate will not harm corneal tissue
- ✓ The high quality, cold-chisel-edges of the carbide steel Alger Brush Burrs are too blunt to harm corneal tissue.
- ✓ Operates with a normal AA battery and low amperage motor means longer battery life
- ✓ Repairable

Codes:

- A4955 Alger Brush Machine
- A4959 Sterile 0.5mm Alger Brush Burrs - Pack of 5
- A4960 Sterile 0.5mm Alger Brush Burrs – Pack of 25



THESE BURRS ARE ALGER Co APPROVED

