




2 Witney Way  
 Boldon Business Park  
 Tyne and Wear  
 NE35 9PE  
 England



# MICROWASH

The following products are now available from  
 Altomed Limited.



To place an order contact Customer Services on  
 Tel: (0) 191 519 0111 Fax: (0) 191 519 0283

**THE Bespoke solution**

Microsurgical instruments require special accommodation during washer machine cleaning, autoclaving and off-site transportation.

MicroWash trays enable effective cleaning and sterilisation whilst also protecting even the finest instruments.

Suitable for all Ophthalmic, Neuro, ENT, Micro-Vascular and Dental instruments. MicroWash trays can be wrapped or placed into rigid containers during transit.

Some examples of tray sizes - Other sizes available

Many carcass sizes are available both in single and double layer configurations.

Sizes that can accommodate a single diamond knife up to a complete DCR, major ear or neurosurgery set give you a wide choice

A9104  
195x70x40mm



A9106  
220x125x40mm



A9110  
250x220x40mm




A9114  
440x250x80mm




**Customising components**



Even the most delicate instruments such as these retinal forceps can be protected during automated washing and off-site processing by accommodating them in a customised MicroWash tray



Each MicroWash tray is supplied with a plastic laminated 'Placement Sheet' which enables both theatre and theatre services staff to relocate instruments into their proper place in the tray



**Scope Trays Also available!**

