



2 Witney Way  
Boldon Business Park  
Tyne and Wear  
NE35 9PE  
England



# FLUORON<sup>®</sup>

## Brilliant Peel Now available in the UK!

The following products are now available from Altomed Limited.

To place an order contact Customer Services on  
Tel: (0) 191 519 0111 Fax: (0) 191 519 0283


**New Higher Density Version Now  
Available...Faster...Higher Intensity**



- Ready to Use
- CE Marked
- Non Toxic
- Stable Stated Concentration

Brilliant Peel selectively stains the Inner Limiting Membrane, and it does so within a short contact time.

The recently added feature of using heavier than standard water in the composition helps overall staining effectiveness



Cat No G-81005

Supplied in a box of 5 individually packed 0.5ml vials

## An I.L.M. selective licensed posterior segment membrane dye to aid peeling