

Tubes can be made to any size required!

Implantable Tubes And Accessories

All products



Marked

NEW

FROSTED FINISH - Performs equally as well as a standard finish tube however it also reduces the possibility of extrusion, the main complication of traditional CDCR.

Due to the multitude of different sizes and styles available, it is not practical to quote every single code number here.

Therefore when ordering a tube please quote the necessary code number by using the format outlined below.


The code is made by quoting first the **Style** (see drawings), the **Flange Width** and the **Length**. If a **Suture Hole** is required add **H** to the end of the code, if it needs to be **Angled** add an **A** and for a **Frosted Finish** add **F**. For example:

Examples - Product Required	Style	Flange Size	Tube Length	With Suture Hole	Angled	Frosted	Code
Lester Jones 3.5mm flange 9mm long with hole and angled	LJ	35	09	H	A		LJ3509HA
Gladstone Putterman 3mm flange 20mm long with hole and straight	GP	30	20	H			GP3020H
Naugle Tube 4mm flange 25mm long no hole and angled	NG	40	25		A		GP3020A
Callahan Cox Tube 3mm flange 28mm long no hole and straight	CC	30	28				CC3028
Lester Jones 3.5mm flange 15mm long, hole, frosted finish	LJ	35	15	H		F	LJ3515HF

Note: All products are supplied as non-sterile as standard

Product Styles And Accessories Available:

[LJ - Lester Jones Tubes, Straight / Angled](#)



2.2mm Outside Diameter

[Jones Gold Dilator](#)



[GP - Gladstone / Putterman Tube](#)



[Polyethylene Tube](#)



Polyethylene Tubes are not Autoclavable

[BY - Baylis 60° Circle Tube](#)



[DC - Double Collar](#)



[CC - Callahan \(Cox\) Tube](#)



3.0mm Outside Diameter

[NG - Naugle Tube](#)



Product Style	Code
Lester Jones Tube	LJ
Lester Jones Tube - Frosted	LJ****F
Callahan Cox Tube (minimum 4.5mm flange)	CC
Gladstone Putterman Tube	GP
Baylis 60° Circle Tube	BY
Naugle Tube	NG
Double Collar Tube	DC
Polyethylene Tube (not autoclavable)	A11908
Accessories - See Also Cox Promotional Sheet	Code
Lester Jones Gold Dilator No. I: 1.7mm and 1.9mm arms	A 11937
Lester Jones Gold Dilator No. II: 2.2mm and 2.4mm arms	A 11937A
Lester Jones Gold Dilator No. III: 2.7mm and 2.9mm arms	A 11937B
Blue anodized aluminium autoclavable sterilisation tray	A 11932
Plastic storage tray with labelled compartments - not autoclavable	A 11930
Pack of 5 Jones nylon cleaning rods - not autoclavable	A 11936
Blue anodized aluminium identification scale for tubes	A 11934
Baylis eyelid plate (globe protector) autoclavable	A 11938